



Coreview

Mammography Phantom

A kit of test objects designed to be used for an ongoing check of imaging performance on dedicated core biopsy imaging systems, and specifically sized to fit the market-leading COREVISION system.

The kit consists of:

- Threshold contrast phantom containing:
 - 6 x Low-contrast, large diameter circular masses
 - 10 x High-contrast, small diameter circular masses
 - 8 x Fibrils
- Spatial resolution test pattern (8-20 LP/mm)
- Geometric distortion phantom
- Uniformity phantom:
 - 1 x 1mm PMMA
 - 1 x 2mm PMMA
 - 2 x 3mm PMMA
- AEC check phantom:
 - 1 x 1mm PMMA
 - 1 x 2mm PMMA
 - 2 x 3mm PMMA




www.leadstestobjects.com



Coreview

Mammography Phantom

Components



- A** Threshold contrast phantom
- B** Uniformity phantom
- C** Spatial resolution test pattern
- D** Geometric distortion phantom