



MagIQ Diffusion

MRI Phantom

The MagIQ Diffusion MRI Phantom is an image quality phantom for calibration of diffusion-weighted (DWI) and apparent diffusion coefficient (ADC) Magnetic Resonance Imaging (MRI) comprising a main chamber and seven 'test tubes' which are supplied filled with different solutions to create a total of 7 separate signal compartments.

Measure:

Contrast values for Diffusion Weighted Imaging (DWI) and Apparent Diffusion Coefficient (ADC) imaging protocols

External dimensions:

220mm diameter
190mm deep

Test Tube dimensions:

35mm diameter

