



TO DR

Dynamic Range Phantom

These test objects are designed for systems which operate by sensing mean signal levels, the TO DR tests dynamic range by measuring detail visibility over a wide range of input Xray intensities. Different designs cater for subtraction and non-subtraction systems.

The TO DR (4 piece)

Supplied in a set that allows it to be used for both subtraction and non-subtraction systems. Comprises of the following test objects:

- Filter plate
- Non-subtraction detail plate
- Subtraction detail plate
- Blank plate

The TO DR (2 piece)

For use with non-subtraction systems only. Comprises of the following test objects:

- Filter plate
- Non-subtraction detail plate



www.leadstestobjects.com